



For a healthy
start in life



12-CHANNEL BODY ARRAY COIL 1.5 / 3.0T

- 12-channel Array, consisting of 2 elements, can be used separately or together as appropriate, comprised of:
 1. part integrated in patient's bed (8-channels)
 2. part as separate coil (4-channels)
- High-resolution MR Imaging of the body of newborns and premature babies, with a weight of up to approx. 4,500g
- Compatible with
 - Siemens (VB, VC, VD, VE)
 - Philips (≥ 3.2 , $\geq 5.1.7$)
 - GE (DV 24/25/26)
 - GE (SV 25.1)

Further compatibilities are in progress

- Approved for use inside MR Diagnostics Incubator System nomag[®] IC *ADVANCED*

Change of technical specification without notice



Kindly provided by:
Clinical Hospital Center, Zagreb, Croatia