



### Incubator with natural air circulation

- Bench top model
- Natural convection
- Entirely made of stainless steel
- Glass door (option Inox door insulated)
- Electronic control
- Door lockup with key
- Heat insulated with fiberglass
- Temperature range from **30°C to 70°C**
- Time regulation from 0 to 600 minutes and infinite time
- Digital display for temperature and time
- Resolution 0,1°C
- Visual and ringing alarms for :
  - Visual and ringing alarms for :
    - over temperature
    - temperature fall
    - probe damage
    - broken heating element
    - black out
- Voltage 230V AC (110V on demand)
- With solid Schuko plug

Art.	Internal dimensions (mm) Dimensions interieur (mm)	External dimensions (mm) Dimensions exterieur (mm)	Capacity (Liters) Capacité (Litres)	Max power (W) Puissance max (W)	Removable shelves Plateaux extractibles
<b>2431/V</b>	385 x 235 x 218 (h)	560 x 370 x 380 (h)	20	100	2



**MeMed**

Advancing health

**MeMed CZ s.r.o.**

Dukelská 5779

43001 Chomutov

Česká republika

Tel. +420 777 333 452

office@memed.cz

www.memed.cz