

Product Data Sheet

Model name: MR5501 (Serie III)

Description: Height adjustable transport bed for the MRI area

7 Tesla tested

Patient weight: up to 220kg

Advantages:

- Hydraulic double pedal pump for adjusting the height
- Side hinged safety railings
- Flame-retardant mattress made of wipeable material and with welded seams
- Adjustable back
- Smooth-running wheels with parking brakes

Dimensions:

- Height: 495 mm to 940 mm
- Mattress thickness: 50 mm

Top:

- Length: 2175 mm
- Width: 640 mm

Bottom:

- Length: 1475 mm
- Width: 610 mm

