

MULTIPAD PRODUCT-FAMILY

THE “NEW STANDARD OF CARE”

Pearltec products represent the new “Standard of Care” when it comes to the placement, positioning and fixation of patients during MR and CT examinations, with regard to image quality, handling, hygiene and patient comfort.

SWISS TECHNOLOGY INSIDE

Together with the ETH Zurich and other renowned institutes, an innovative and patented product technology has been developed consisting of soft and mobile polystyrene beads combined with inflatable air chambers, which fits the patient perfectly, thereby creating individual and gentle placement, positioning and fixation.

INNOVATIVE SOLUTIONS

Pearltec offers its customers innovative solutions for placement, positioning and fixation in radiology and radiotherapy. This includes a wide range of standard products, tailor-made customer solutions for systems suppliers and technical solutions for research.



MULTIPAD FAMILY OF PRODUCTS


The flexible and universally applicable MULTIPAD family of products enables optimum patient placement, positioning and fixation when it comes to examinations using head and extremity coils.



MULTIPAD Standard

Thanks to its dimensions of 17x6 cm and the variable thickness of 2 to 5 cm, the MULTIPAD Standard is particularly suitable for situations where there are no specific restrictions.

MR



MULTIPAD Plus

Thanks to its dimensions of 17x11 cm and the variable thickness of 2 to 5 cm, the MULTIPAD Plus is particularly suitable to fill larger gaps.

MR



MULTIPAD Slim

Thanks to its dimensions of 17x9 cm and the variable thickness of 1 to 3 cm, the MULTIPAD Slim is particularly suitable for situations where space is limited.

MR



MULTIPAD cover


The covers for the MULTIPAD product family form an integral component for a hygienic and patientfriendly application.

MR

At least four good reasons to consider Pearltec products as the new “Standard of Care” for the placement, positioning and the fixation of your patients for MRI and CT examinations.


FEWER MOTION ARTEFACTS

Leads to improved image quality and therefore a reliable diagnosis. Leads to fewer repetitions and therefore saves time and costs.



USER-FRIENDLY SYSTEM SOLUTION

Enables efficient, clear and straightforward handling, therefore reducing the large variety of different placement products.




PERFECTLY HYGIENIC

It is easy to clean, provides reliable protection and therefore fosters the safety and wellbeing of the patient.



HIGH LEVEL OF PATIENT COMFORT

The pressure-free, adaptive immobilisation increases patient comfort and therefore has a positive effect on the entire examination procedure.



APPLICATION FIELD

The flexible MULTIPAD family, consisting of MULTIPAD Standard, MULTIPAD Plus, MULTIPAD Slim can be universally used either alone or in combination during all examinations with head and extremity coils.

