


# Revolution Front Elastic Back 403



- Cool Wear moisture-wicking interior fabric
- Professionally designed for an exceptional fit and comfort
- SmartTrack enabled



Model Number	Description
403	Revolution Front Elastic Back
Colors	500, 600, 700, 800, 900, Remix
Core	Lightweight Lead KIARMOR Bi-Layer Low Lead KIARMOR Bi-Layer Lead-Free
Protection	.25mm LE, .35 mm LE, .50mm LE
Sizes	Male S-5XL Female XS-5XL Length alterations ±2" or ±4"
Fabric	Revolution Fabric with Cool Wear
Optional Base	B303, B404 Classic Nylon Colors without Cool Wear (Custom Alterations)
Monogram	<b>Block</b> or <i>Script</i>

### Radiation Protection Options



STANDARD CONFIGURATION		OPTION 1		OPTION 2	
.50mm LE		.35mm LE		.25mm LE	
Absorbed Dose	Attenuation	Absorbed Dose	Attenuation	Absorbed Dose	Attenuation
0.9%	99.1%	2.45%	97.6%	4.9%	95.1%
0.6%	99.4%	1.9%	98.1%	3.5%	96.5%
0.6%	99.4%	2.0%	98.0%	3.9%	96.1%
Highest Protection		Best Combination of Protection & Weight		Lightest Weight	

Shown in 615 Black/Orange

Model	Description	Gender	Sizes	Color	Protection Options		
					Lightweight Lead	KIARMOR Low Lead	KIARMOR Lead-Free
403	Revolution Front Elastic Back	M	S, M, L, XL, 2XL, 3XL, 4XL, 5XL	Revolution	LL	KI-LL	KI-LF
		F	XS, M, L, XL, 2XL, 3XL, 4XL, 5XL	Revolution	LL	KI-LL	KI-LF