

**Vest & Skirt**  
**VAS**



- Full overlapping front panels for maximum protection
- Hook and loop closures
- Double welt pockets
- Comfortable shoulder pads
- Belt with buckle closure on skirts
- Skirt hanging straps



Shown in SH-4 Stylized Hearts

Model Number	Description
VAS	Vest & Skirt
Colors	Classic Nylon Colors
Core	Lightweight Lead KIARMOR Bi-Layer Low Lead KIARMOR Bi-Layer Lead-Free
Protection	.50mm LE at overlap on front, .25mm LE on back
Sizes	Male S-5XL Female XS-5XL
Monogram	<b>Block</b> or <i>Script</i>

Model	Description	Gender	Sizes	Color	Protection Options		
					Lightweight Lead	KIARMOR Low Lead	KIARMOR Lead-Free
VAS	Vest & Skirt	M	S, M, L, XL, 2XL, 3XL, 4XL, 5XL	Classic Nylon Colors	LL	KI-LL	KI-LF
		F	XS, S, M, L, XL, 2XL, 3XL, 4XL, 5XL	Classic Nylon Colors	LL	KI-LL	KI-LF