


Thyroid Collar without Binding

NST



Shown with Magnetic Closure Option (-M)

Dimensions | Hook and Loop & Magnetic Closures



Product	Thyroid Collar without Binding
Features	Comfortable fit Large protection area Optional monogram
Sizes	X-Small, Standard, Medium
Colors	Revolution 700 series & All Classic Nylon Colors
Core	Lightweight Lead KIARMOR Bi-Layer Low Lead KIARMOR Bi-Layer Lead-Free
Protection	.25mm LE, .35 mm LE, .50mm LE
Closure	Hook and Loop Closure, Buckle, Magnet
Fabric	Revolution or Classic
Monogram	Block or <i>Script</i>

Model Number	Description
NST-XS	Thyroid Collar without Binding, Size X-Small, Hook and Loop Closure
NST-XS-M	Thyroid Collar without Binding, Size X-Small, Magnetic Closure
NST-XS-B	Thyroid Collar without Binding, Size X-Small, Buckle Closure
NST-S1	Thyroid Collar without Binding, Size Standard, Hook and Loop Closure
NST-S1-M	Thyroid Collar without Binding, Size Standard, Magnetic Closure
NST-S1-B	Thyroid Collar without Binding, Size Standard, Buckle Closure
NST-M1	Thyroid Collar without Binding, Size Medium, Hook and Loop Closure
NST-M1-M	Thyroid Collar without Binding, Size Medium, Magnetic Closure
NST-M1-B	Thyroid Collar without Binding, Size Medium, Buckle Closure