



Wall Mounted Apron Rack
683429



- Heavy gauge steel
- 200 lbs. Capacity
- Powder coated
- Scratch resistant
- Easy to install
- Requires only minimal space
- Stores up to 10 aprons
- Optional vest and skirt adapter (Part Number 683000)

Model Number	Description
683429	Wall Mounted Apron Rack
Dimensions	12" H x 36" W X 26" D
Colors	White, Custom Options
Options	Powder Coated Steel Stainless Steel (Upgrade)
Shipping	Box 1: 40 lbs (16" x 58" x 9")