



For a healthy
start in life





8-CHANNEL HEAD ARRAY COIL 1.5 / 3.0T

- 8-channel Array, cylindrical, receive only coil
- High-resolution MR Imaging of the brains of newborns and premature babies with a head circumference of up to approx. 40cm
- Compatible with
 - Siemens (VB, VC, VD, VE)
 - Philips (≥ 3.2 , $\geq 5.1.7$)
 - GE (DV 24/25/26)
 - GE (SV 25.1)

Further compatibilities are in progress

- Approved for use inside MR Diagnostics Incubator System nomag[®] IC *ADVANCED*

Change of technical specification without notice



Kindly provided by:
Oslo University Hospital, Rikshospitalet, Norway