

# Revolution Full Wrap Black Belt

## 203



Shown in 616 Black/Volt



Hook and loop closure side slits

- Partial overlap to increase protection
- Hook and loop closure side slits for increased mobility
- Cool Wear moisture-wicking interior fabric
- Professionally designed for an exceptional fit and comfort
- Snap on, fully adjustable Black Belt®
- SmartTrack enabled

Model Number	Description
203	Revolution Full Wrap Black Belt
Colors	500, 600, 700, 800, 900, Remix
Core	Lightweight Lead KIARMOR Bi-Layer Low Lead KIARMOR Bi-Layer Lead-Free
Protection	.25mm LE, .35mm LE, .50mm LE
Sizes	Male S-5XL Female XS-5XL Length alterations ±2" or ±4"
Fabric	Revolution Fabric with Cool Wear
Optional Base	B203 Classic Nylon Colors without Cool Wear (Custom alterations)
Monogram	<b>Block</b> or <i>Script</i>

### Radiation Protection Options

STANDARD CONFIGURATION		OPTION 1		OPTION 2	
.50mm LE		.35mm LE		.25mm LE	
Absorbed Dose	Attenuation	Absorbed Dose	Attenuation	Absorbed Dose	Attenuation
0.9%	99.1%	2.45%	97.6%	4.9%	95.1%
0.6%	99.4%	1.9%	98.1%	3.5%	96.5%
0.6%	99.4%	2.0%	98.0%	3.9%	96.1%
Highest Protection		Best Combination of Protection & Weight		Lightest Weight	



Back panel is always .25mm LE unless requested otherwise

Model	Description	Gender	Sizes	Color	Protection Options		
					Lightweight Lead	KIARMOR Low Lead	KIARMOR Lead-Free
203	Revolution Full Wrap Black Belt	M	S, M, L, XL, 2XL, 3XL, 4XL, 5XL	Revolution	LL	KI-LL	KI-LF
		F	XS, M, L, XL, 2XL, 3XL, 4XL, 5XL	Revolution	LL	KI-LL	KI-LF