



Magnetic Thyroid Collar

MTC



- Comfortable fit
- Large coverage area
- Optional monogram
- Optional KTC attachment

Dimensions | Hook and Loop Closures & Magnetic Closures



Model Number	Description
MTC	Magnetic Thyroid Collar
Size	X-Small, Standard, Medium
Colors	All
Core	Lightweight Lead KIARMOR Bi-Layer Low Lead KIARMOR Bi-Layer Lead-Free
Protection	.25mm LE, .35 mm LE, .50mm LE
Closure	Magnet
Fabric	Revolution or Classic
Monogram	Block or <i>Script</i>

Model Number	Description
MTC-XS	Magnetic Thyroid Collar, Size X-Small, Magnetic Closure
MTC-S1	Magnetic Thyroid Collar, Size Standard, Magnetic Closure
MTC-M1	Magnetic Thyroid Collar, Size Medium, Magnetic Closure